# **Transform Your Clinic with Skin Analyser 3D**

Elevate your consultations and improve patient outcomes with the advanced Skin Analyser 3D. Its cutting-edge technology ensures accurate, engaging diagnostics, helping you deliver tailored treatments with confidence

# SKIN ANALYSER 3D



Learn more about Aesthetic Consults information @aesthetic\_consults official account

## Contact us

Shop 23A Mount Annan Market place 11-13 main st, Mount Annan 2567

+61 489 904212; +61 2 8526 4838









At **Aesthetics Consults**, we bring you the cutting-edge **Skin Analyser 3D**, a revolutionary diagnostic device that uses **3D imaging and AI-powered technology** to provide in-depth insights into your clients' skin. With **8-spectrum imaging** and advanced analysis tools, this device enables clinics to deliver highly customised treatments with precision and confidence

### Why Choose Skin Analyser 3D for Your Clinic?

Skin Analyser 3D transforms the way you approach skin diagnostics, making your consultations more accurate, engaging, and professional

### Key Features –

**3D Imaging Technology:** Creates a detailed 3D model of the face for accurate assessments

8-Spectrum Imaging: Detects a wide range of skin concerns, including pigmentation, acne, and sensitivity

**AI-Powered Analysis:** Rapidly generates precise results, tailored to each client's unique skin profile

**Comprehensive Reporting:** Provides easy-to-read reports for tracking progress over time

**User-Friendly Design:** Intuitive touchscreen interface with quick image capture and analysis

Learn more about Aesthetic Consults information @aesthetic\_consults official account





# **How Skin Analyser 3D Works**

Using a combination of advanced imaging technology and AI algorithms, the Skin Analyser 3D evaluates a range of skin concerns in real time, providing actionable insights for treatment planning



#### Controlled lock



Three images



Comprehensive report



✓ Single detection



✓ 8 Screen Display



✓ A single report



- 3D Face Mapping: Creates an accurate, high-resolution 3D model of the face for visual analysis
- 8-Spectrum Imaging: Utilises different light spectrums to detect surface and subsurface concerns
- Al-Driven Diagnosis: Processes detailed data to identify
- Customised Reports: Generates detailed client reports with easy-to-understand visuals and personalised treatment recommendations

### **Comprehensive Diagnostic Capabilities**

8-Spectrum Imaging: The Skin Analyser 3D uses advanced light spectrums to assess various skin concerns

- White Light: Analyses surface texture and pore size
- UV Light: Highlights clogged pores and acne
- Polarised Light: Maps shallow and deep pigmentation
- Brown Light: Detects age spots and melasma
- Red Light: Evaluates sensitivity and inflammation
- Mixed Light: Measures hydration and moisture levels

Real-Time 3D Visualisation: Engage clients with detailed 3D models of their face, allowing them to visualise their concerns and understand the recommended treatments

Progress Tracking: Compare results over time with easy-to-read reports that document improvements in wrinkles, pigmentation, and hydration

Progress Tracking: Compare results over time with easy-to-read reports that document improvements in wrinkles, pigmentation, and hydration

issues such as wrinkles, pigmentation, pores, and hydration

• Oil-Rich Light: Assesses sebum levels for oily or acne-prone skin

# Benefits of Adding Skin Analyser 3D to Your Clinic

### **For Your Clinic**

- Improve Client Confidence: Build trust with data-driven, visual evidence of skin concerns
- Expand Service Offerings: Combine diagnostics with facials, laser treatments, or skincare plans
- Boost Retention: Progress tracking encourages clients to stay committed to their treatment plans
- Maximise Revenue: Offer premium diagnostic services and increase client satisfaction

## **For Your Patients**

- Customised Treatment Plans: Tailored recommendations based on detailed skin analysis
- Engaging Visuals: Clear and detailed 3D imaging helps clients understand their skin better
- Trackable Results: Tangible proof of improvement with each follow-up session
- Safe and Non-Invasive: Suitable for all skin types, including sensitive skin



### **Device Specifications**

- Camera Resolution: 20-megapixel high-definition imaging
- 8 Imaging Spectrums: Detects a wide range of concerns with precision
- Screen: Intuitive 15.6-inch colour touchscreen for easy operation

- reports

# **Frequently Asked Questions (FAQ)**

A full session takes approximately 15-20 minutes, including imaging and report generation

### Can it track progress over multiple treatments?

Yes, the device provides comparative reports to monitor improvements in real time

### Is it suitable for all skin types?

Absolutely! The Skin Analyser 3D is safe and effective for all skin tones and types

### Do I need special training to operate it?

The device is user-friendly, and we provide comprehensive training to ensure your team is confident using it

### Can the analysis be combined with other treatments?

Yes, it pairs perfectly with facials, chemical peels, laser treatments, and other advanced skincare options

• Compact Design: Lightweight and ergonomic (400 x 300 x 350mm, 10kg) Connectivity: USB and Bluetooth-enabled for seamless integration • Report Generation: Compatible with A4 colour printers for professional